

DATA SHEET LATEX GLOVES MAXLATICECS 10 PF

DESCRIPTION AND FEATURES



LEGAL REGULATIONS

Class I: Medical Device according to Regulation (EU) 2017/745, with standards EN 455-1, EN 455-2, EN 455-3. EN 455-4.

Category III: Personal Protective Equipment (PPE) according to Regulation (EU) 2016/425, with standards EN 21420:2020, EN 374-1:2016+A1:2018 Type B, EN 374-4:2013, EN 374-5:2016.

MaxLaticecs 10 PF disposable gloves are made of natural latex without powder, specially designed to provide maximum protection in demanding environments. Their superior thickness ensures exceptional puncture resistance, making them an ideal choice for tasks requiring increased protection. Suitable for medical, therapeutic, and industrial applications.

- Made from natural latex without powder.
- Micro-textured surface at the fingertips.
- High resistance for a high level of safety.
- Ambidextrous.
- AQL 1.5 for consistent quality.
- Shelf life: 5 years.

INSTRUCTIONS FOR USE

It is important to note that latex gloves can cause allergic reactions. If you experience any allergies, consult a doctor immediately.

Before putting on the gloves, make sure they have no defects; insert your hand inside and adjust them carefully. Store them in a dry and cool place, away from light.

They are intended for single use and should be disposed of in an environmentally friendly manner in appropriate waste containers, in accordance with current regulations.

FEATURES

Deference	C:	Davi	Box Carron	Length		Minimum thickness in mm	
Reference	Size	вох		in mm		Finger	Palm
37002	S	100 units	10 boxes	240	80 ±10	0,22 ± 0,02	0,20 ± 0,02
37003	М	100 units	10 boxes	240	95 ± 10	0,22 ± 0,02	0,20 ± 0,02
37004	L	100 units	10 boxes	240	110 ± 10	0,22 ± 0,02	0,20 ± 0,02
37005	XL	100 units	10 boxes	240	≥ 110	0,22 ± 0,02	0,20 ± 0,02

Before aging

	Resistence	Elasticity			
	Min 18 Mpa.	650%			
Después del envejecimiento					
	Resistence	Elasticity			
	Min 14 Mpa.	500%			

CONTACT US

Visit our website



Calle los Verdiales 28, Valdemoro (Madrid) (+34) 910 376 792 - administracion@icemedics.com