



Technical Information:

Disposable visitor kit for minimal risk environments and single-use applications.

- Included Products
 - Ref. 5300B: White Polyethylene LDPE Gown with Snap Closures
 - Ref. 2100B: White Accordion-Style Polypropylene Cap
 - Ref. 6300A: Extra-Rough Blue CPE Polyethylene Shoe Cover
 - Ref. 3500A: Blue Type IIR Surgical Mask with Ear Loops
- SINGLE USE ONLY
- Latex-Free Product

Standards and Usage:

- Complies with EU Regulation 2017/745 on medical devices Class I
- Complies with EU Regulation 2016/425 on personal protective equipment PPE Category I
- Complies with **EU Regulation 10/2011** on plastic materials intended for food contact, in accordance with Article 5, Paragraph 1 of Regulation (EC) No. 1935/2004
- The mask meets the requirements of EN 14683:2019+AC2019 standards.
- Manufactured under ISO 13485:2003+AC2007 standards

Usage Instructions:

- 1. Verify the product's condition and expiry date before use.
- 2. Wash hands thoroughly before handling.
- 3. Unpack the visitor kit and carefully wear the components, ensuring proper adjustment in the designated areas for optimal protection.
- 4. Do not reuse.

Logistics Information:

Labeling Compliant with Applicable Standards

Presentation:

- Packaging: Bag containing 1 unit with compliant labeling
- Master Carton: 60 units
- Carton Dimensions: 35 x 29 x 30.5 cm
- Euro Pallet: 48 cartons

| Ref. | Color | EAN Code |
|--------|------------|---------------|
| IC-KIT | Blue&White | 0685071749364 |

Store in a dry, clean place at temperatures between 5°C and 35°C. Avoid direct exposure to sunlight.

