

TECHNICAL DATA SHEET

Laticecs PF



PRODUCT INFORMATION

- ✓ Made of powder-free natural latex.
- ✓ Micro-textured surface on fingertips for better grip.
- ✓ High strength for a high level of safety.
- ✓ Ambidextrous for convenient use by right and left-handed users.
- ✓ AQL 1.5 for consistent quality.
- ✓ Lifetime: 5 years



Medical



Food industry



Industry



Automotive



Cleaning


APPLICATIONS AND USES

- ✓ It is important to note that latex gloves can cause allergic reactions. In case of allergy, consult a doctor immediately.
- ✓ Before putting on the gloves, make sure they are free of defects; place your hand inside and gently pull them on.
- ✓ Store them in a cool, dry place, protected from light.
- ✓ They are intended for single use only and must be disposed of in an environmentally friendly manner in appropriate waste disposal containers in accordance with current regulations.

LEGAL REGULATIONS

Class I: Medical Device according to EU Regulation 2017/725, whose standards are: EN 455-1, EN 455-2, EN 455-3, EN 455-4 and ASTM D-6978-5.

Category III: Personal Protective Equipment (PPE) according to EU Regulation 2016/425, whose standards are: EN 21240:2020, EN 374-1:2016+A1: 2018 Type B, EN 374-4:2019, EN 374-5:2016.



TECHNICAL SPECIFICATIONS

| Reference | Size | Length in mm | Width in mm | Minimum thickness in mm | |
|-----------|------|--------------|-------------|-------------------------|-------------|
| | | | | Finger | Palm |
| 28012 | S | 240 | 85 ± 5 | 0,12 ± 0,02 | 0,10 ± 0,02 |
| 28013 | M | 240 | 95 ± 5 | 0,12 ± 0,02 | 0,10 ± 0,02 |
| 28014 | L | 240 | 105 ± 5 | 0,12 ± 0,02 | 0,10 ± 0,02 |
| 28015 | XL | 240 | 110 ± 5 | 0,12 ± 0,02 | 0,10 ± 0,02 |

| Before aging | |
|--------------|------------|
| Resistance | Elasticity |
| 18 Mpa. | 650% |
| After aging | |
| Resistance | Elasticity |
| 18 Mpa. | 500% |

TECHNICAL DATA SHEET

Laticecs PF



LOGISTICS



100 gloves
per box



10 boxes
per carton



70 cartons
per pallet

TARIFF CODE: 4015.19

MINIMUM ORDER QUANTITY: CARTON OF 10 DISPENSERS PER BOX OF 100 GLOVES

| Unit of consumption: dispenser or box of 100 units | | | | |
|--|-----------------------------|----------------------------|--------|-------|
| Code and description | Weight box / dispenser (Gr) | Box / dispenser dimensions | | |
| | | Long | Width | High |
| 28012- Laticecs PF S.S | ± 540 Gr | 215 mm | 120 mm | 70 mm |
| 28013- Laticecs PF S.M | | | | |
| 28014- Laticecs PF S.L | | | | |
| 28015- Laticecs PF S.XL | | | | |

| Unit of sale: 10 dispensers or box per carton | | | | |
|---|--------------------|----------------------------|--------|--------|
| Code and description | Carton weight (Kg) | Box / dispenser dimensions | | |
| | | Long | Width | High |
| 28012- Laticecs PF S.S | ± 5,40 Kg | 360 mm | 245 mm | 225 mm |
| 28013- Laticecs PF S.M | | | | |
| 28014- Laticecs PF S.L | | | | |
| 28015- Laticecs PF S.XL | | | | |

| Transport unit: Pallet | | | | | | |
|------------------------|---------------------------|----------------|--------------------|------------------|--------|---------|
| Number of cartons | Number of cartons per row | Number of rows | Pallet weight (kg) | Dimension Pallet | | |
| | | | | Long | Width | High |
| 70 cartons | 10 cartons | 7 rows | ± 415 kg | 1200 mm | 800 mm | 1700 mm |

| Description | EAN-13 Inner | GTIN-14 Outer |
|--------------------------------|---------------|----------------|
| Laticecs PF, Latex glove S. S | 8437020302074 | 18437020302071 |
| Laticecs PF, Latex glove S. M | 8437020302081 | 18437020302088 |
| Laticecs PF, Latex glove S. L | 8437020302098 | 18437020302095 |
| Laticecs PF, Latex glove S. XL | 8437020302104 | 18437020302101 |